

AHN myLab® SM-01

Single Plate Magnetic Stirrer 1500 rpm

Usage - in laboratories, pharmacies and schools, mixing liquid industry

Maximum speed - from 15 to 1500 rpm with an interval of 50 rpm / Software driven gradual acceleration

Time setting - ranging from 1- 99 minutes / continuous operation feature

Improved mixing - pulse mode where the mixing direction changes every 30 seconds

Maintenance free design - water spill, corrosion and chemically resistant / includes no moving parts / high quality stirring bar

User-friendly interface and design - digital display showing the real-time speed and the remaining time / last run memory function / extremely thin (only 16.0 mm thick)/ flat surface enables easy wiping and cleaning



Ordering information

Description

AHN myLab® SM-01 Single Plate Magnetic Stirrer 1500 rpm

Cat. No.

7-300-00-0

Delivery package

- 1 pc. Magnetic Stirrer 1500 rpm
- 1 pc. Stirrer bar (PTFE coated, 25 x 8 mm)
- 1 pc. Instruction manual
- 1 pc. Warranty card

Specifications

Max. stirring quantity	800 mL
Recommended stir bar length	25 mm
Protection class	IP 65
Dimensions (WxDxH)	120 x 164 x 16 mm
Weight	400 g



+49(0)3631/46594-04
+49(0)3631/46594-10
info@ahn-bio.de

Uthleber Weg 14
D-99734 Nordhausen
www.ahn-bio.de

ahn®
biotechnologie GmbH