

AHN myLab® VT-03

Digital Vortex Mixer 4200 rpm

Usage - comfortable everyday mixing operations

Maximum Speed - adjustable up to 4200 rpm

Digital display - adjustable speed setting with an interval of 10 rpm and time setting from 0 - 999 minutes / displays speed and time alternatively at every 5 seconds

Stable operation - through inbuilt counter balance system

Robust design - brushless DC motor enables maintenance free and long lasting operation


Programmable pulse mode - ON time and OFF time can be individually modified

Multiple applications - optional attachments for micro plate, micro tubes and flask (***only applicable in combination with the universal platform 7-110-01-0**) / sticker for permissible speed for various attachments

Small footprint saves valuable bench space - suitable for use in hoods and cold rooms



Optional attachments *



Ordering information

Description

AHN myLab® VT-03 Digital Vortex Mixer 4200 rpm

Cat. No.

7-120-00-0

Delivery package

- 1 pc. Digital Vortex Mixer 4200 rpm
- 1 pc. Power supply adaptor
- 1 pc. Standard cup head (pre-installed)
- 1 pc. Instruction manual
- 1 pc. Warranty card

Optional attachments

Cat. No.

- Universal platform **7-110-01-0**
- Foam for microtubes* **7-110-02-0**
- Foam for microplate* **7-110-03-0**
- Rubber holder for flask* **7-110-04-0**

*only applicable in combination with the universal platform 7-110-01-0.

Specifications

| | |
|-----------------------|----------------------|
| Orbital diameter | 4 mm |
| Different modes | Touch, continuous |
| Speed setting | 300-4200 rpm |
| Max. loading capacity | 500 g |
| Dimensions (WxDxH) | 205 x 136 x 138.5 mm |
| Weight | 2.8 kg |

