

AHN myLab® VT-02

Vortex Mixer 2800 rpm

Usage - everyday mixing operations

Maximum Speed - adjustable up to 2800 rpm

Comfortable use - selection of different modes (touch, stand-by and continuous), depending on the desired use / smooth and silent operation

Stable operation - through inbuilt counter balance system

Multiple applications - optional attachments for micro plate, micro tubes and flask (***only applicable in combination with the universal platform 7-110-01-0**) / sticker for permissible speed for various attachments

Small footprint saves valuable bench space - suitable for use in hoods in the temperature range from 5 to 40 °C



Optional attachments *



Ordering information

Description

AHN myLab® VT-02 Vortex Mixer 2800 rpm

Cat. No.

7-110-00-0

Delivery package

- 1 pc. Vortex Mixer 2800 rpm
- 1 pc. Power supply adaptor
- 1 pc. Standard cup head (pre-installed)
- 1 pc. Instruction manual
- 1 pc. Warranty card

Optional attachments

Cat. No.

- | | |
|--------------------------|------------|
| Universal platform | 7-110-01-0 |
| Foam for microtubes* | 7-110-02-0 |
| Foam for microplate* | 7-110-03-0 |
| Rubber holder for flask* | 7-110-04-0 |

Specifications

Orbital diameter	4 mm
Speed accuracy	± 5.0 %
Max. loading capacity	500 g
Protection class	IP 21
Power supply	220 - 240VAC, 50Hz
Dimensions (WxDxH)	205 x 136 x 138.5 mm
Weight	3.0 kg

*only applicable in combination with the universal platform 7-110-01-0.



+49(0)3631/46594-04
+49(0)3631/46594-10
info@ahn-bio.de

Uthleber Weg 14
D-99734 Nordhausen
www.ahn-bio.de

ahn®
biotechnologie GmbH